

zeen Adjustment Tables

Lift Ranges in Pounds (kg)						
LIFT SETTINGS	Medium Springs		Lighter Springs		Heavier Springs	
1	100 - 130	(45 - 59)	50 - 75	(23 - 34)	150 - 185	(68 - 84)
2	130 - 160	(59 - 73)	75 - 100	(34 - 45)	185 - 220	(84 - 100)
3	160 - 190	(73 - 86)	100 - 125	(45 - 57)	220 - 255	(100 - 116)
4	190 - 220	(86 - 100)	125 - 150	(57 - 68)	255 - 290	(116 - 131)
5	220 - 250	(100 - 113)	150 - 175	(68 - 79)	290 - 325	(131 - 148)

Standard Stand-Over Height Ranges by Setting in Inches (cm)		
Seated	21"	(53)
1	23" - 24"	(58 - 61)
2	24" - 25"	(61 - 64)
3	25" - 26"	(64 - 66)
4	26" - 27"	(66 - 69)
5	27" - 28"	(69 - 71)
6	28" - 29"	(71 - 74)
7	29" - 30"	(74 - 76)
8	30" - 31"	(76 - 79)
9	31" - 32"	(79 - 81)
10	32" - 33"	(81 - 84)
11	33" - 34"	(84 - 86)
12	34" - 35"	(86 - 89)
13	35" - 36"	(89 - 91)

Tall Extension Settings - To be purchased for those that exceed the standard maximum stand-over height - All settings are 2" higher than the equivalent setting on the standard chart		
Seated	23"	(58)
1	25" - 26"	(64 - 66)
2	26" - 27"	(66 - 69)
3	27" - 28"	(69 - 71)
4	28" - 29"	(71 - 74)
5	29" - 30"	(74 - 76)
6	30" - 31"	(76 - 79)
7	31" - 32"	(79 - 81)
8	32" - 33"	(81 - 84)
9	33" - 34"	(84 - 86)
10	34" - 35"	(86 - 89)
11	35" - 36"	(89 - 91)
12	36" - 37"	(91 - 93)
13	37" - 38"	(93 - 96)



NOTE: Zeen Seat Cushion adds an additional inch of support under the user.



Lift Adjustment



Height Adjustment



Settings Review

